


NO	DESC	PICTURE
Aksesoris Bronchoscopy		
1	Mucus Aspirator Medbar , Penampung cairan bilasan, 1 pc	
2	MPC-ST, Mouthpiece Standart Fuji , Reusable, Box/5	
3	MPC-P, Mouthpiece Pediatric Fuji , Reusable, Box/5	
4	Mouthpiece Standart with strap Medbar , Disposable, Box/10	
5	BF1812SF, Biopsy Forcep Broncho Standart Fenestrated Fuji , Working Channel 2.0 mm, Effective Length 1200mm. Reusable, 1 pc	
6	GF1812V, V-Shape Grasping Forcep Bronchoscopy Fuji , Working Channel 2.0 mm, Effective Length 1200mm. Reusable, 1 pc	
7	BIO1-C4-18-120, Biopsy Forceps Bronchoscopy Medwork , Oval Jaws, Coated, Diameter Jaws 1.8mm, Working Channel 2.0 mm, Effective Length 120cm. Disposable, 10 pc/Box	
8	WB3212FW2, Channel Washing Brush for Scope Broncho Fuji , 1 bristle , Working Channel 2.0 ~ 2.8mm, Effective Length 1200mm. 1 pc	


Aksesoris Endoscopy (Gastro & Colon)		
9	BIO1-C4-18-180, Biopsy Forceps for Gastro Medwork , Oval Jaws, Coated, Pediatric, Diameter Jaws 1.8mm, Working Channel 2.0 mm, Effective Length 180cm. Disposable, 10 pc/Box	
10	BIO1-C4-23-180, Biopsy Forceps for Gastro Medwork , Oval Jaws, Coated, Diameter Jaws 2.3mm, Working Channel 2.8 mm, Effective Length 180cm. Disposable, 10 pc/Box	
11	BIO1-C4-23-230, Biopsy Forceps for Colon Medwork , Oval Jaws, Coated, Diameter Jaws 2.3mm, Working Channel 2.8 mm, Effective Length 230cm. Disposable, 10 pc/Box	
12	BIO2-C1-23-180, Biopsy Forceps for Gastro Medwork , Oval Jaws, Diameter Jaws 2.3mm, Working Channel 2.8 mm, Effective Length 180cm. Reusable, 1 pc.	
13	BIO2-C1-23-230, Biopsy Forceps for Colon Medwork , Oval Jaws, Diameter Jaws 2.3mm, Working Channel 2.8 mm, Effective Length 230cm. Reusable, 1 pc.	
14	BF2416SF, Biopsy Forceps for Gastro Standart Fenestrated Fuji . Working Channel 2.8 mm, Effective Length 1600mm. Reusable, 1 pc	
15	BF2424SF, Biopsy Forceps for Colon Standart Fenestrated Fuji . Working Channel 2.8 mm, Effective Length 2400mm. Reusable, 1 pc	
16	GF2416A, Grasping Forcep Alligator Jaw For Gastro Fuji , Working Channel 2.8 mm, Effective Length 1800mm. Reusable, 1 pc	
17	GF2424A, Grasping Forcep Alligator Jaw For Colon Fuji , Working Channel 2.8 mm, Effective Length 2400mm. Reusable, 1 pc	

18	SX-H-2300, Shangxian Hemoclip Handle (Delivery System) , Sychro rotation, Working Channel 2.8mm, Working Length 2300mm. Reusable, 1 pc	 
19	SX-C-135, Hemoclip , Jaw Angle 135°, working length 5mm, max opening 12mm, closing length 13.5mm. Disposable. 30 pc/ Box.	 
20	SX-C-90, Hemoclip , Jaw Angle 90°, working length 5.7mm, max opening 12mm, closing length 14.2mm. Disposable. 30 pc/ Box.	 
21	WB3218FW2, Channel Washing Brush for Scope Gastro Fuji , 1 bristle , Working Channel 2.0 ~ 2.8mm, Effective Length 1800mm. 1 pc	
22	WB3521FW2, Channel Washing Brush for Scope Colon Fuji , 1 bristle , Working Channel 2.8 ~ 3.2mm, Effective Length 2100mm. 1 pc	

Jarum Core Biopsy

1	RBG-100 Möller Medical Blue Biopsy Gun , Reusable, automatic biopsy instrument for use with disposable biopsy needles to obtain soft tissue samples	
2	DNG-1020 Disposable Biopsy Needle , compatible with Möller Medical Blue and Promag™ Ultra	
3	COX - Coaxial Needle , to guide the biopsy needles and to allow multiple sampling with a single puncture	
4	ABG-2020 Automatic Biopsy Gun , Automatic disposable biopsy device for obtaining high quality histological samples	

5	SBG-2010 Semi-automatic Biopsy Gun , Semi-automatic disposable biopsy device for obtaining high quality histological samples	
6	Chiba Needles FNB-1001 , Fine needle biopsy, For percutaneous extraction of cytological and histological samples	
7	Spinal Needles FNB-1002 , Fine needle biopsy, For percutaneous extraction of cytological and histological samples	
8	BLM-1110 Breast Lesion Marker , Repositionable Breast Lesion Marker for tagging and pre-operative localisation	
9	BLM-1120 Stiff Marker and BLM-1130 Flexible Marker , Breast Lesion Marker for tagging and pre-operative localisation	
10	BLM-1125 Stiff Marker and echogenic tip , BLM-1135 Flexible Marker and echogenic tip , Breast Lesion Marker for tagging and pre-operative localisation	
11	BOW-1060 Bone Marrow Biopsy System , Needle set for atraumatic bone marrow biopsy	

12	BOW-1070 Bone Marrow Aspiration Biopsy , Needle set for atraumatic bone marrow biopsy and aspiration	 The image shows a BOW-1070 Bone Marrow Aspiration Biopsy needle set. It consists of a blue handle with a white grip, a long needle, and two red caps.
13	BOY-1080 Bone Biopsy Needles , Needle set for atraumatic bone biopsy	 The image shows a BOY-1080 Bone Biopsy Needles needle set. It consists of a blue handle with a white grip, a long needle, and two red caps.

FUJIFILM

medwork

 **MöllerMedical**

medbar®


SHANGXIAN