

「Cleaner and Healthier Air in the World」

Plasma Air Sterilizer

DWK-DP3000



Distributor
INNOSONIAN


Dawoo Korea
Manufacturer

Air needs cleaning

As you wash dirty clothes, the air you breathe needs to be cleaned.

All-in-one purification, sterilization, disinfection, and deodorization of dirty indoor air with fine particles, bacteria, virus, odor, and other elements that contribute to sick building syndrome!

Dawoo Korea will lead a clean Korea with Plasma Air Sterilizer, which is embedded with innovation including air pressure plasma, photocatalytic sterilization and disinfection technologies, an antibacterial HEPA filter, and carbon absorption technologies.



Plasma Air Sterilizer

DWK-DP3000



Pre-vacuum Air Removal

Touch Screen type Steam Sterilizer

- ▶ **A New Selection for Large Steam Sterilizer**
- ▶ **Realized the Ergonomic Design for Maximizing Operator's Safety and Convenience.**



Automatic Vertical Sliding Door

HS-5035PT/HS-5025PT/HS-5020PT/HS-4085PT

DOUBLE DOOR type



SINGLE DOOR type



SINGLE & DOUBLE Door type



HS-5035PT



HS-5025PT



HS-5020PT



HS-4085PT

It is very convenient to operate and control the sterilizer by light touching the graphic designed button easily and is able to monitor the cycle progress in real time on the large color touch screen.

The door moves up and down smoothly without any physical force by adopting the automatic vertical sliding system, therefore the large size sterilization is much more safe and convenient than the hinged door type. Moreover it has a advantage of saving the required space for installing and operating the sterilizer.